



Power up your IPR

Precision Swiss-made, diamond-impregnated files accomplish a wide range of tasks including breaking contact, reduction, contouring, fine trimming and finishing.

SpaceFile® Interproximal Reduction Files

- Can be used three ways: Hand Held, with Removable Grip, or in optional Air-driven Handpiece
- Diamond-Impregnated, double-sided or single-sided
- Designed to allow low temperature operation and prevent clogging
- Files, grip and inner tray are autoclavable
- Color coded to match Endodontic file system

Contra-Angle Handpiece

- Use with all SpaceFile® Files
- Compact design makes reaching posterior regions easy
- Use with most standard slow-speed air motors & compressed air units such as A-dec, KaVo, W&H, NSK, and Midwest
- Reciprocating motion of 1.6mm
- Autoclavable
- One-Year Warranty

Assorted Packs

Includes 8 Files (1 of each thickness), 1 Serrated Contact Point Saw, and 1 Removable Grip

Double-Sided

Single-Sided
Left

Single-Sided
Right

40034

40034-LS

40034-RS

Single Type Packs

Includes 8 files (all same thickness) and 1 Removable Grip

THICKNESS

Polishing

A final step in tooth surface restoration.



40035-G

40035-LS-G

40035-RS-G

0.085mm

0.070mm

0.070mm

Finishing

For finishing prior to final polish.



40035-P

40035-LS-P

40035-RS-P

0.110mm

0.100mm

0.100mm

XX Fine

Pre-polish and finish reduced areas. Also use as initial strip for excessive crowding.



40035-W

40035-LS-W

40035-RS-W

0.140mm

0.123mm

0.123mm

X Fine

Pre-polish stripped areas. Also use as initial strip for crowded cases.



40035-Y

40035-LS-Y

40035-RS-Y

0.160mm

0.138mm

0.138mm

Fine

Use to initiate reduction in crowded cases. Also use to contour after reduction.



40035-R

40035-LS-R

40035-RS-R

0.270mm

0.180mm

0.180mm

Medium

Use for reduction with some crowding or when some space exists. Smaller grain can achieve results before using Coarse strip.



40035-BL

40035-LS-BL

40035-RS-BL

0.330mm

0.210mm

0.210mm

Coarse

Primary working strip for IPR. Also use to trim and contour.



40035-GR

40035-LS-GR

40035-RS-GR

0.370mm

0.235mm

0.235mm

X Coarse

Use for maximum reduction at start of IPR.



40035-BLK

40035-LS-BLK

40035-RS-BLK

0.490mm

0.290mm

0.290mm

Serrated Contact Point Saw

Use to break contact points on enamel, composite, amalgam, porcelain or gold.



40035-CS

—

—

SpaceFile® Removable Grip

Pack of 3

40035-GRIP

Contra-Angle Handpiece (sold separately)

40027ANG

Air Motor (sold separately)

40027MOT



A-dec, KaVo, W&H and NSK are not trademarks of Dentsply Sirona. Colors are for illustration purposes only. Actual color of product may vary.

SpaceFile®: Invented by Dr. Carlos Navarro. Design protected by: US Patent 0057540A1; EU Patent 4341095. MediMark® Europe BP 2332, F-38033 Grenoble Cedex 2, France. Rx Only. Made in Switzerland.

Dentsply Sirona Orthodontics Inc. • 7290 26th Court East • Sarasota, FL 34243 USA

40032-MKT Rev. 06 05-2021